



Marianas LightSheet

Versatile Multiview Light Sheet Microscopy
System for Imaging Model Organisms



Marianas LightSheet

Marianas LightSheet (MLS) combines the low phototoxicity and large sample space of dual inverted selective plane illumination microscopy (diSPIM) with the power and flexibility of a research-grade inverted microscope system. MLS supports a wide variety of specimen types and sizes as well as experimental designs including photoablation, photomanipulation, and computer-generated holography. Add spinning disk confocal, TIRF and FLIM to transform the system into a versatile multimodal imaging system all controlled with SlideBook software.



TTL Sync


Millisecond timing and trigger
Control of multiple devices


Scan-Optimized X,Y Stage

High-precision movement of sample during long-term experiments

LaserStack

Modular laser combiner
Up to eight lasers






Fiber Switcher

Up to four fiber outputs
Millisecond path switching

Additional Light Path

Add cameras, photomanipulation devices,
spinning disk confocal and more



Spinning Disk Confocal

Live cell 3D confocal imaging
Super-resolution dual microlens disk




Ablate!

Laser ablation system
355nm and 532nm


Vector2

Epi mounted or camera port
mounted photomanipulation
Modular high-speed X,Y scanner



Vector3

Epi-mounted TIRF and photomanipulation
Spinning X,Y TIRF with expansive FN22 FOV
Liquid light guide input for LED light source



Light Sheet Scanners

High-speed MEMs mirror or cylindrical
lens scanners

Environmental Control

Stage top and cage incubators
Temperature, gas and humidity



Motorized SPIM Z Drive

High precision multi-position and vertical montage capture

sCMOS Cameras

95% QE low read noise
Large field-of-view options

Objective Piezos


150µm or 300µm piezos for long-term high-speed 3D sampling

Objectives

Water and multi-immersion dual-view optimized objective pairs


Fully Automated Research Microscope

Motorized objective, condenser and path selection
Autofocus (Definite Focus 3)
PSF-optimized objectives




Scanning Modes

Fixed
Ultra-fast single plane time-lapse capture with stationary light sheet.




ZEBRAFISH HEART
Zebrafish heartbeat (GFP-labeled Cardiomyocytes). Courtesy of Dr. Jamie Nichols, University of Colorado.

Piezo
Fast 3D capture up to 3 volumes/sec of smaller specimens with synchronized sheet and objective movement.




ARABIDOPSIS ROOT
Piezo slice scan capture of an Arabidopsis root (YFP-labeled membrane). Courtesy of Dr. David Ehrhardt, Carnegie Institution for Science.

Stage
Large 3D capture up to 1 volume/sec of specimens translated through a stationary light sheet.




CLEARED BRAIN SLICE
Stage scan captures of a cleared mouse brain (GCaMP-GFP) with a pair of 10x (top) and 40x (bottom) objectives. Courtesy of Dr. Rob Campbell, University College London.

Dual-View Piezo
Alternating Piezo captures combined with SlideBook Multiview Reconstruction for isotropic resolution.



DICTYOSTELIUM
Dual-view piezo slice scan timelapse of Dictyostelium cells (membrane-rhodamin).

Dual-View Stage
Alternating Stage captures combined with SlideBook Multiview Reconstruction for isotropic resolution.




ZEBRAFISH
Dual-view stage scan timelapse (GFP-actin). Courtesy of Elric Esposito, Spencer Shorte UtechS Photonic Biolmaging and Nicolas Dray, Laure Bally-Cuif Zebrafish Neurogenetics lab, Institut Pasteur.

Objective Pairs


Magnification	Numerical Aperture	Refractive Index	Maximum Field of View	Recommended Use
Water Immersion				
10x	0.3	1.33-1.4	1490µm	Largest field of view
20x	0.5	1.33-1.4	620µm	Balanced field of view and resolution
40x	0.8	1.33-1.4	370µm	Highest resolution
Multi-Immersion				
10x	0.3	1.33-1.56	1470µm	Largest field of view
16x	0.4	1.33-1.56	850µm	Higher resolution large field
24x	0.7	1.33-1.56	590µm	Highest resolution

Flexible Sample Mounting



Dual Inverted Selective Plane Microscopy




Dual inverted selective plane microscopy (diSPIM) uses dual camera and scanner objective pairs positioned at 45° to the specimen plane to alternate capture and excitation between **Path 1** and **Path 2**. The flexible system geometry allows for single- and dual-sided imaging with conventional sample mounting. SlideBook seamlessly controls all hardware modalities, deskewing, stitching/montaging, Multiview Reconstruction and joint deconvolution resulting in isotropic sub-cellular resolution across a wide range of samples. Marianas LightSheet integrates diSPIM with a research-grade inverted microscope enabling complex multi-modal light sheet experiments with photoablation, photostimulation and computer-generated holography.



ZEBRAFISH
Dual-view stage scan capture of a Zebrafish (GFP-Actin). Path 1 (left) and path 2 (right). Courtesy of Dr. Cody Smith, University of Notre Dame.

SlideBook Multiview Reconstruction


SlideBook for MLS supports several methods to deskew, rotate, register and deconvolve diSPIM images. Dual-view capture registration is performed either through Point of Interest Detection or 3i's proprietary Local Cross-Correlation and Registration algorithms. Results are instantly displayed in 3D for quick adjustment to achieve an ideal fit before joint deconvolution produces a fused image. Large datasets that would otherwise require high-powered computer resources can be automatically split into a computer's available RAM and computed in parallel utilizing GPU-optimized processes.



Photomanipulation

Ablate! Laser Ablation System


- 355nm or 532nm pulsed laser
- Fixed point or galvo-scanned variable region of interest
- 2D or 3D regions
- Diffraction limited spot



PHOTOABLATION
Mammalian cell with membrane labeled in green. Before (left) and after ablation (right).

Vector2 Scanning Photomanipulation


- Photoactivation/FRAP
- Galvo-scanned variable region of interest
- 2D or 3D regions
- Diffraction limited spot



PHOTOBLEACHING
Mammalian cell with membrane labeled in green before bleaching and FRAP curve in SlideBook.


Phasor Holographic Photomanipulation

- Spatial light modulator-generated holography for optogenetics stimulation/FRAP/voltage imaging
- Simultaneous 3D stimulation of multiple, separate regions
- Visible and multiphoton stimulation without scanning




3D PHOTOSTIMULATION
3D illumination pattern (left) applied to a 3D specimen (right) to stimulate multiple regions simultaneously.

Total Internal Reflection Fluorescence



Vector3 Epi-Mounted TIRF and Photomanipulation

Vector3 is a motorized spinning X,Y TIRF system with built-in photomanipulation and widefield fluorescence capabilities. Vector3 offers an expansive TIRF field-of-view (FN22) designed for modern sCMOS cameras. Spinning the TIRF beam via galvo mirrors results in smooth, evenly illuminated images across the field. Vector3 is capable of galvo-scanned photomanipulation of diffraction-limited spots via user-drawn ROIs for easy FRAP and photoconversion experiments. Incorporation into the epi path of a microscope allows the use of any existing microscope cameras – for example, cameras attached to a CSU-W1 spinning disk – keeping the alignment identical between different imaging modalities and enabling powerful, multimodal imaging experiments. Integrated collimating optics are compatible with most LED illuminators (via a liquid light guide) for widefield illumination.



Data courtesy of Dr. Chris Bakal and Oliver Inge at the Institute of Cancer Research in London.

Spinning Disk Confocal

Yokogawa spinning disk confocals utilize a dual Nipkow disk with microlenses for the best optical sectioning and minimal pinhole crosstalk. This proven technology is the best solution for intravital imaging where optical sectioning and speed are both critically important.



CSU-X1

- Highest speed imaging at up to 2000fps
- Field of view 7mm x 10mm
- 50µm pinhole disk with microlenses
- Manual and motorized versions





CSU-W1

- High speed imaging up to 200fps
- Wide field of view 16mm x 17mm
- 25µm and 50µm pinhole disks for lower and higher magnification objectives
- Motorization including disk exchange, variable aperture, camera port selection and camera port magnification
- Options for split-view imaging, NIR imaging, illumination field flattening and super-resolution imaging

CSU-W1 SoRa

CSU-W1 SoRa is an easy-to-use super-resolution microscopy solution utilizing a dual Nipkow disk assembly with microlenses on both the illuminating and pinhole disks. SoRa images have a 1.4x resolution improvement and deconvolved SoRa images have a 2x resolution improvement compared to standard spinning disk data. With a maximum speed of 200fps, low phototoxicity and no limitation on dyes or fluors, SoRa is ideal for super-resolution intravital imaging. SoRa is also available as an upgrade to existing CSU-W1 systems.




Imaging of microtubules in fixed bovine pulmonary artery endothelial cells. Azuma, T. and Kei, T. (2015) Super-resolution spinning-disk confocal microscopy using optical photon reassignment. Opt Express. Jun 1;23(11):15003-11. doi: 10.1364/OE.23.015003.

Uniformizer | Illumination Field Uniformity


For exceptionally even illumination across the entire field, Uniformizer conditions the gaussian beam from the illumination fiber optic to distribute light evenly across the field. Using a set of microlens arrays, Uniformizer flattens the field to as little as 1% variance and boosts overall intensity up to 50%.



Standard Montage



With Uniformizer





SlideBook software supports research microscopy through the entire experimental process. By managing everything from instrument control to image processing and data analysis, SlideBook allows scientists to focus on investigation rather than instrumentation. SlideBook controls hundreds of instruments in and around the microscope from dozens of manufacturers enabling researchers to integrate their preferred components and upgrade to the latest devices once available.

User-Selectable App Appearance

Select a color scheme from dozens of options
Switch on-the-fly from dark to light themes

SlideBook Open File Format

Directory-based open file format for big data and high performance computing applications

Multiview Reconstruction

SlideBook can use a number of methods to deskew, rotate, register and deconvolve images captured with MLS. Data is visualized, deskewed and rotated with a single click. Initial shifts and registration results can be instantly checked in 3D and adjusted before joint deconvolution to a fused image.

Volume Rendering

3D and 4D visualization tools rotate and deskew MLS data on-the-fly and support a user specified bounding box and storyboard interface where multiple perspectives can be assembled into a single movie.

NVIDIA CUDA GPU Acceleration



GPU acceleration of computationally-intensive operations such as deconvolution

Multi-Position and Montage

Capture single timepoints or time-lapses at multiple positions with different experimental setups and throughout the range of capture modes

3D Capture Status

Volumetric projection during 4D capture supported across all instruments



Capabilities



Capture

Control hundreds of devices including microscopes, stages, lasers, wheels, piezos, scanners, shutters and much more.



View

Visualize data through any numbers of portals, from single images to z-stacks, time lapse, color channels and 4D views.



Analyze

Analyze images and extract statistical data via a wide variety of algorithms while maintaining original data integrity.



Scripting

Macro scripting for capture and analysis enhances the flexibility and power available to users.



Communicate

Present and export data easily as 16-bit TIFFs, 3D movies, graphs or spreadsheets. Data is directly portable to MATLAB and Excel and adheres to Open Microscopy Environment (OME) standards.

Partners



MATLAB

Through hierarchical and conditional capture, user-supplied MATLAB programs can control experimental workflows.



Aivia

Aivia is an innovative and complete 2D-to-5D image visualization, analysis and interpretation platform with artificial intelligence-guided image analysis.



syGlass

syGlass enables 3D and 4D visualization and analysis of SlideBook data in a virtual reality environment.



Dell

The latest high-power computer workstations control all microscope hardware and enable high-speed processing, segmentation and volume rendering of terabyte (TB) datasets.

Systems Engineering

3i's Systems Engineering department designs, builds and extensively tests every customer system. From spinning disk confocal to multiphoton to lightsheet and photomanipulation, 3i has delivered over a thousand custom, cutting-edge microscopy systems to help answer some of the most complex scientific questions.



Application Knowledge | Scientific Consulting

A team of PhD scientists meet with each client to document and better understand the scientific context of the user group to ensure that the capabilities of the delivered system match the underlying research goals.

Performance Criteria | Targeted to Experiments

Understanding key experiments and imaging paradigms allows Systems Engineering to apply targeted testing criteria to every system.

Customized Hardware | Novel Light Creation

No matter how complex or customized a light path may be for imaging or photostimulation, our engineers ensure that light is manipulated and directed to where it is needed, when it is needed.

Custom Test Plan | Assure Experiment Success


When a technically advanced experiment requires specific system performance to succeed in the lab, a custom test plan assures the system meets that mark prior to delivery.

System Integration | Synchronization of Dozens of Instruments

Systems Engineering combines institutional knowledge and scientific consultation to ensure that the instruments in each system are configured for experimental success in the lab.

System Test Report | Guaranteed Performance

Performance metrics and results of the custom test plan are documented in a System Test Report delivered with each system.



System Test Report


General System Testing	Initials
Confirm presence of all system components from quoted order	
Verify correct configuration hardware properties, optical paths, and objectives in SlideBook settings	
Confirm parfocality across objectives	
Test field flatness	
Generate Point Spread Function for high NA objectives (83x)	

LaserStack	Initials
Verify correct lasers and correct power rating are present	
Verify laser power out of head matches Manufacturing Test Sheet	
Measure laser power at fiber input and output	
Measure laser power at back aperture	
Measure laser power at objective output (all wavelengths with 63x)	
Calculate system throughput	


Motorized XY Stage	Initials
Level stage insert	
Validate motorized travel set distances with stage micrometer	
Double-click to center calibration	
Confirm accuracy of XY multipoint acquisition	

Zeiss Z-Drive	Initials
Test z-position repeatability	


Page 6




System Test Report




Arabidopsis Embryo



Zebrafish Eye




Drosophila melanogaster embryo




Zebrafish Embryo

Page 14


Application Data





DROSOPHILA EYE
Stage scan Y-montage of a Drosophila pupal retina showing utrophin labeled F-actin in pseudocolour. Courtesy of Courtney Lancaster, Pichaud lab, University College London Laboratory for Molecular Cell Biology (LMCB).




CYANOBACTERIA
Maximum intensity projection (MIP) of a stage scan timelapse capture of cyanobacteria biofilm autofluorescence excited with 488nm and 561nm after 48 hours. Courtesy of the Institute of Cancer Research, London.




BIOFILM
Dual-view stage scan of Candida albicans growth (GFP-Mitochondria). Courtesy of Dr. D Pentland, collaboration Gourley-& Laissue-Labs, Universities of Kent & Essex.



ZEBRAFISH TAIL
Dual-view stage scan capture of a Zebrafish (GFP-Actin). Courtesy of Dr. Cody Smith, University of Notre Dame.





HUMAN EMBRYONIC KIDNEY (HEK) CELLS
Dual-view piezo scan capture of HEK cells with labeled tubulin (green), actin (magenta) and nuclei (blue). Courtesy of Dr. Deirdre Kavanagh, University of Birmingham.




DROSOPHILA
Dual-view stage scan capture of a Drosophila wing disc showing Viking (magenta) and membrane (green). Courtesy of Dr. Yanlan Mao and Dr. Rob Tedley, University College London Laboratory for Molecular Cell Biology (LMCB).


Application Data




 **MOUSE BRAIN**
Stage scan capture of iDISCO-cleared embryonic mouse brain with antibody stainings for smooth muscle cells (white) and endothelial cells (magenta). Courtesy of Tijana Perovice, Gerhardt lab, Max-Delbrück Center for Molecular Medicine.





ARABIDOPSIS EMBRYO
Dual-view piezo scan capture of an Arabidopsis embryo showing autofluorescence excited with 405nm (blue) and 488nm (green). Courtesy of University of Warwick.





ZEBRAFISH EMBRYO
Stage scan capture of a zebrafish embryo 20 hours post-fertilization with labeled actin, tubulin and nuclei.




 **DROSOPHILA EMBRYO**
Piezo scan timelapse of synchronised cell division of Drosophila melanogaster embryo (~ 4hours old, H2B-labeled chromatin).















 **ARABIDOPSIS TRICHOMES**
Piezo scan timelapse of Arabidopsis trichomes (GFP labeled actin).



 **DROSOPHILA EMBRYO**
Stage scan capture of Drosophila embryo with labeled nuclei (orange) and corpora allata cells (blue). Courtesy of Dr. Lucas Dent, the Institute of Cancer Research, London.

Support and Maintenance

A variety of software and equipment support levels help keep systems running well for years. A Software Support Agreement allows labs to run the latest version of SlideBook with new acquisition and analysis features. It includes direct access to 3i staff via email, phone and video chat. A System Maintenance Agreement adds an annual preventative maintenance visit, 3i service visits and 3i coordination of any repairs, although repair and replacement parts are not included. A System Extended Warranty adds full coverage for repairs and replacement parts. Additionally, 3i application scientists may provide in-person and webinar-based application training.

	Software Maintenance	System Maintenance	System Warranty
Phone, Email and Video Chat Support			
SlideBook Software Releases			
Service Visits and Annual PM Visit			
Repairs Coordinated by 3i			
Application Training In-Person or Online			
Full Warranty Coverage of all System Hardware			

BUILT BY SCIENTISTS FOR SCIENTISTS

3i designs and manufactures technologies for living cell, live cell, and intravital fluorescence microscopy including superresolution, computer-generated holography, spinning disk confocal, multi-photon and lightsheet. SlideBook software manages everything from instrument control to image capture, processing and data analysis. 3i was established in 1995 by a group of cell biologists, neuroscientists, and computer scientists to provide advanced multi-dimensional microscopy platforms that are intuitive to use, modular in design, and meet the evolving needs of investigators in the biological research community.



USA

Intelligent Imaging Innovations, Inc.
3575 Ringsby Court, Suite 102
Denver, CO 80216 USA
+1 303 607 9429

UK

Intelligent Imaging Innovations Ltd
17 Westbourne Studios
242 Acklam Rd, London, UK W10 5JJ
+44 2037 440031

DE

Intelligent Imaging Innovations GmbH
Königsallee 9-21
D-37081 Göttingen, Germany
+49 551 508 39 266

FR

Intelligent Imaging Innovations SAS
17 rue Moreau
75012 Paris, France